

Water Washable Dental Resin

Suitable for dental restoration, invisible correction thermoforming and other applications.

- Washable
- Low Skin Irritation
- Low Skin Allergenic
- Low Shrinkage
- Temperature Resistance
- High Hardness



Learn More

Color: **Dark wheat skin**
 Specification: **500/1000g**
 Compatible with **DLP/LCD**
 Applicable light source band: **385-405nm**



*1.Quality Management System (ISO - 13485) in the process of certification

Resin Property (Post-cure)

Hardness	ASTM D 2240	85-90	Shore D
Flexural modulus	ASTM D 790	1000-1500	MPa
Flexural strength	ASTM D 790	60-75	MPa
Tensile modulus	ASTM D 638	600-1000	MPa
Tensile strength	ASTM D 638	40-55	MPa
Elongation at break	ASTM D 638	10-15	%
Impact strength,notched Izod	ASTM D 256	35-45	J/m
Heat deflection temperature	ASTM D648 @66PSI	75-80	°C