


Removable Die Model Resin

Suitable for the fabrication of dental working and situation models, control models and orthodontic models.

- Low skin irritation
- Low skin allergenic
- Low shrinkage
- Temperature resistance
- High hardness physical properties



Learn More

Color: **Wheat Complexion**
 Specification: **500/1000g**
 Compatible with **DLP/LCD**
 Applicable light source band: **385-405nm**



*1.Quality Management System (ISO - 13485) in the process of certification

Resin Property (Post-cure)

Hardness	ASTM D 2240	85-90	Shore D
Flexural modulus	ASTM D 790	1500-2000	MPa
Flexural strength	ASTM D 790	60-80	MPa
Tensile modulus	ASTM D 638	600-1200	MPa
Tensile strength	ASTM D 638	35-60	MPa
Elongation at break	ASTM D 638	10-15	%
Impact strength,notched Izod	ASTM D 256	35-50	J/m
Heat deflection temperature	ASTM D648 @66PSI	75-85	°C