

Surgical Guide Resin Pro

Suitable for Dental and Orthopedics surgery fixture and medical application.

- Low Skin Irritation
- Low Skin Allergenic
- Low Shrinkage
- High Transparency
- High Elongation
- High Impact Strengyh



Learn More

Color: **Transparent**
 Specification: **500/1000g**
 Compatible with **DLP/LCD**
 Applicable light source band: **385-405nm**



Resin Property (Post-cure)

Hardness	ASTM D 2240	75-80	Shore D
Flexural modulus	ASTM D 790	1000-1200	MPa
Flexural strength	ASTM D 790	>40	MPa
Tensile modulus	ASTM D 638	230-270	MPa
Tensile strength	ASTM D 638	>30	MPa
Elongation at break	ASTM D 638	110-140	%
Impact strength,notched Izod	ASTM D 256	240-300	J/m
Heat deflection temperature	ASTM D648 @66PSI	60-70	°C

Withstands high-temperature sterilisation at 135°C without cracking or deformation.